

EASE ADJUST BY THERAFIRM®


30-50
mmHg*



INELASTIC COMPRESSION BELOW KNEE WRAP

Adjustable gradient inelastic compression wrap for venous and lymphatic disease.

- ✓ **ANATOMICAL DESIGN** Contours to the calf for a better fit, plus seamless straps reduce pressure points
- ✓ **EASY APPLICATION** Molded hook for easy peel straps, and an inner support sleeve secures wrap in position while donning and adjusting
- ✓ **DAY OR NIGHT USE** Inelastic materials offer a high working pressure during the day and a low resting pressure with strong containment at night
- ✓ **MAXIMUM DURABILITY** Made with 3mm thick premium grade materials to offer strong containment for managing venous edema or lymphedema
- ✓ **LINER SOCKS INCLUDED** A pair of 10-15 mmHg Core-Spun® by Therafirm Socks provide compression and added comfort in the foot and ankle
- ✓ **PRODUCT GUARANTEE** Therafirm provides a 6-month warranty for materials and workmanship



SIZE	ANKLE	CALF	UNDER KNEE
Small	6.5-9.5" (16-24cm)	9.5-14" (24-35cm)	10-14" (25-36cm)
Medium	8-11" (20-28cm)	12.5-17" (32-43cm)	12.5-17" (32-43cm)
Large	9.5-12.5" (24-32cm)	16-20" (40-51cm)	15.5-19.5" (39-50cm)
X-Large	11-14" (28-36cm)	19-23" (48-59cm)	18-22.5" (46-57cm)
REGULAR LENGTH Less than 14" (36cm)		TALL LENGTH More than 14" (36cm)	

Ankle, calf, and length are required for proper fit. Sizes based on average.

*Designed to provide 30-50 mmHg pressure based on wearer needs.

CONTRAINDICATIONS: Compression products should not be worn and are contraindicated if you have any of the following conditions: arterial insufficiency or degeneration, acute DVT (deep vein thrombosis), untreated congestive heart failure, untreated cancer, untreated infection, absent or severely impaired sensation, allergy to materials, severe or moderate peripheral arterial disease, severe cognitive impairment.

wholesale.therafirm.com



CONNECT WITH US!

5P-EA-ADJS 2306

THERAFIRM®
THUASNE

HOW TO MEASURE + SIZE



MEASURING FOR THE EASE ADJUST WRAP

1. On the affected leg, measure the circumference of the patient's ankle at the narrowest point, directly above the ankle bone.
2. Measure the circumference of the calf at the widest point.
3. Measure the circumference under the patient's knee, just below the knee bend.
4. Measure the length from the floor to the knee bend (in a straight line).

Use the table below to find the size that corresponds to the measurements. If any measurement is between sizes, select the larger size.

Ankle, calf, and length are required for proper fit. Sizes based on average.

TIPS FOR GETTING THE PERFECT FIT

Patients should be measured while standing when they first get up in the morning or after elevating their legs and feet for 15 minutes. This will help prevent measuring based on the swelling of the affected area. Fitters should take snug measurements against the skin while patients are standing. Fitters should never remove any wound dressings and avoid measuring over socks or slacks.

MEASUREMENTS	SMALL		MEDIUM		LARGE		X-LARGE	
	REGULAR	TALL	REGULAR	TALL	REGULAR	TALL	REGULAR	TALL
C UNDER KNEE	68615	68715	68625	68725	68635	68735	68645	68745
B CALF	10–14" (25–36cm)		16–20" (40–51cm)		15.5–19.5" (39–50cm)		18–22.5" (46–57cm)	
A ANKLE	9.5–14" (24–35cm)		12.5–17" (32–43cm)		16–20" (40–51cm)		19–23" (48–59cm)	
D LENGTH	6.5–9.5" (16–24cm)		8–11" (20–28cm)		9.5–12.5" (24–32cm)		11–14" (28–36cm)	
	Regular Less than 14" (36cm) Tall More than 14" (36cm)		Ankle, calf, and length are required for proper fit. Sizes based on average.					

HOW TO ORDER
Phone 800-821-3094
Email cs@knitrite.com



CURRENT CUSTOMERS
Scan to order online via your wholesale portal account.



NEW CUSTOMERS
Scan to fill out our wholesale inquiry form.